

**1. Aim:**

Transiently transfect Expres2 cells in 24-well format for production of secreted recombinant proteins

**2. Applicability:**

Sequence-verified expression constructs in Expres2 vector for expression in S2 cells

**3. Responsibilities:**

Main Responsible: Christian Aigner

Substitute: Orla Dunne

**4. Required materials:**

Chemically competent E. coli NEB10 $\beta$  cells (self-made)

LB, SOC media (media kitchen)

EX-CELL<sup>®</sup> 420 Serum-Free Medium for Insect Cells (Sigma 14420C)

24 well block (Qiagen 19583)

DNA purification mini-prep kit (AnalytikJena 845-KS-1041250)

0.2  $\mu$ M syringe filter (VWR 514-0073)

Corning 250 ml Flask vented cap (Corning 431406)

Expres2 Insect – transfection reagent TRx5 (Expres2ionbio S2-55A-001)

500mM CuSO<sub>4</sub>

**5. Procedure:**

**A. DNA Preparation**

Prior to transfection, sufficient amounts of plasmid DNA need to be prepared (7.5  $\mu$ g), this can be prepared in advance and stored at -20° C.

**Step 1**

Transform plasmids into NEB10 $\beta$  cells (200  $\mu$ l of cells split into 4) using the standard procedure. Reconstitute cells in 300  $\mu$ l of SOC media for 1 h at 37°C and split into 2 x 5 ml LB+Antibiotic (mini A and mini B) and place at 37° C overnight (hands on time for 24 transformations: 1.5 h)

**Step 2**

Mini prep one 5 ml culture per plasmid (mini A) using the standard procedure, however, for the last wash and elution steps do this using the bunsen burner (2.5 h).

**Step3**

Check the concentration, and if required mini-prep mini B. Mini A and B can be pooled.

## B. Transfection

### Day 0 (preferably Wednesday, 15min)

Split the cells by centrifugation (450 g, 3 min.) and resuspend the cells to  $8 \times 10^6$  cells/ml in EX-CELL<sup>®</sup> 420 Serum-Free Medium for Insect Cells in shake flasks and incubate at 25° C at 110 rpm.

### Day 1 (Thursday, 45min)

- Count the cell stock (usually 1:3 dilution necessary)
- Prepare ~80ml with  $8 \times 10^6$  cells/ml in EX-CELL<sup>®</sup> 420 Serum-Free Medium (72ml needed) – centrifuge the cells (450 g, 3 min) and resuspend in 80 ml media.
- Transfer 3ml of the cells for each transfection into a well of the 24-well plate (use Multipette – 10 ml Tip)
- Add 7.5 µg DNA to each transfection
- Add 37.5 µl Express2 Insect – TR 5x to each transfection (use Multipette – 0.5 ml Tip)
- Gently shake the plate to mix the transfection Reaction
- Incubate 5 minutes at room temperature
- Incubate the 24-well plate at 200 rpm at 25 °C for 24 h

### Day 2 (Friday, 5min)

Induce each well with 3 µl 500mM CuSO<sub>4</sub> and incubate at 25° C at 110 rpm for 3 days.

### Day 5 (Monday, 1h)

- Harvest the cultures by centrifugation 20 min at 4000 g at 4 °C.
- Using a syringe with a needle, remove the supernatant from the 24-well plate and filter (0.2 µM syringe filter) it into a 15ml Falcon tube which you can store at 4 °C.